

DNA Markers & Ladders From ID Labs™!

Ladder	Fragment Sizes	Cat#	Size (Gels)
100 bp Ladder, Original	100, 200, 300, 400, 500, 600, 700, 800, 900, 1000,1500	IDBB2593-1	100
100bp-1kb Ladder, Original	100, 200, 300, 400, 500, 600, 800, 1000,1500	IDBB1007-1 IDBB1007-2	100 500
100bp-1kb Ladder, Ready to Use	100, 200, 300, 400, 500, 600, 800, 1000,1500	IDBB1008-1 IDBB1008-2	100 500
DNA Logical Ladder, Original	100, 200, 300, 400,600, 800,1000,1500, 2kb, 3kb, 4kb, 6kb, 8kb,10kb	IDBB1009-1 IDBB1009-2	100 500
DNA Logical Ladder, Ready to Use	100, 200, 300, 400, 600, 800,1000,1500, 2kb, 3kb,4kb, 6kb, 8kb,10kb	IDBB1010-1 IDBB1010-2	100 500
DNA KB Ladder, Original	0.5kb, 1kb, 2kb, 3kb, 4kb, 5kb, 6kb, 8kb,10kb	IDBB1011-1 IDBB1011-2	100 500
DNA KB Ladder, Ready to Use	0.5kb, 1kb, 2kb, 3kb, 4kb, 5kb ,6kb, 8kb,10kb	IDBB1012-1 IDBB1012-2	100 500
Lambda DNA/ Eco RI, Original	3530, 4878, 5643, 5804, 7421, 21226	IDBB2001-1 IDBB2001-2	50 250
Lambda DNA / Hind III , Ready to Use	125, 564, 2027, 2322, 4361, 6557, 9416, 23130,	IDBB2002-1 IDBB2002-2	100 500
Lambda DNA / Hind III PLUS, Ready to Use	125, 564, 330, 1000, 1500, 2027, 2322, 4361, 6557, 9416, 23130	IDBB2003-1 IDBB2003-2	100 500

02/07

Molecular Biology Research Products From ID Labs™!

Product	Cat#	Size	Product	Cat#	Size
X-Gal (Ultra Pure)	ID0083100	100 mg	Proteinase K	IDBB2008-1	100 mg
	ID0083001	1 g		IDBB2008-2	500 mg
	ID0083010	10 g		IDBB2008-3	1 g
IPTG (Ultra Pure)	ID0168001	5 g	dNTP Mix (10mM)	IDBB1001A	0.5 ml
	ID0168005	10 g		IDBB1001B	4 x 0.5 ml
	ID0168025	25 g		IDBB1001C	10 x 0.5 ml
Lambda DNA (dam-, dcm-)	IDBB2004-1	0.1 mg	dNTP Mix (25mM)	IDBB1002A	0.5 ml
	IDBB2004-2	0.5 mg		IDBB1002B	4 x 0.5 ml
Agarose A (Biotech Grade)	IDBB2592-1	500g	d-ATP (100mM)	IDBB1002C	10 x 0.5 ml
	IDBB2592-2	1 kg		IDBB1003A	0.5 ml
Other Agarose Types	Inquire	Inquire	IDBB1003B	4 x 0.5 ml	
RNase A	IDBB2006-1	250 mg	IDBB1003C	10 x 0.5 ml	
(DNase Free, High Purity)	IDBB2006-2	500 mg	d-GTP (100mM)	IDBB1004A	0.5 ml
	IDBB2006-3	1 g		IDBB1004B	4 x 0.5 ml
RNase Inhibitor (Biotech Grade)	IDBB2007-1	2000 U		IDBB1004C	10 x 0.5 ml
RNase A Solution (10mg/ml)	IDBB2007-2	5 x 2000U	d-CTP (100mM)	IDBB1005A	0.5ml
	IDBB2633	1ml		IDBB1005B	4 x 0.5 ml
X-Gluc (Ultra Pure)	IDBB2588-1	100mg		IDBB1006C	10 x 0.5 ml
	IDBB2588-2	500mg	d-TTP (100mM)	IDBB1006A	0.5 ml
		IDBB1006B		4 x 0.5 ml	
		IDBB1006C		10 x 0.5 ml	

Other Molecular Biology Products Available:

- IDPure™ DNA & RNA Purification Kits
- Custom Oligonucleotide Synthesis

Please Email us to Place an Order or inquire about Bulk Prices and More Details!

www.idlabs.com

info@idlabs.com

02/07