



Product Catalogue

CYTOGENETICS _ IDETECTTM PROBES &
TISSUE CULTURE 2010

Our Mission: To provide Life Science, Medical Drug Discovery Specialists with high quality,
innovative and cost-effective immune and gene/chromosome-based medical products and services.



Cytogenetics

Accessories

CounterSafe™ Spray Disinfectant	IS1131	250 ml	Nontoxic, Biodegradable Solution for disinfecting surfaces of lab benches & incubators, Cytogenetics, Accessories	<input type="checkbox"/>
HydroSafe™-I Disinfectant	IS1110	100 ml	Solution for disinfecting of Water Baths in CO2 Incubators (Suggest 10 ml/L), Cytogenetics, Accessories	<input type="checkbox"/>
HydroSafe™-II Disinfectant	IS1111	50 ml	Solution for disinfecting standard water baths (Suggest 2ml/L), Cytogenetics, Accessories	<input type="checkbox"/>
IDBlock-It™ DNA (Bovine)	IDBK1003	0.5 mg	Bovine repeat sequense DNA, Cytogenetics, Accessories	<input type="checkbox"/>
IDBlock-It™ DNA (Human)	IDBK1001	0.5 mg	Human repeat sequense DNA, Cytogenetics, Accessories	<input type="checkbox"/>
IDBlock-It™ DNA (Mouse)	IDBK1002	0.5 mg	Mouse repeat sequense DNA, Cytogenetics, Accessories	<input type="checkbox"/>
IDBlock-It™ DNA (Ovine)	IDBK1006	0.5 mg	Ovine/Sheep repeat sequense DNA, Cytogenetics, Accessories	<input type="checkbox"/>
IDBlock-It™ DNA (Porcine)	IDBK1005	0.5 mg	Porcine/Pig repeat sequense DNA, Cytogenetics, Accessories	<input type="checkbox"/>
IDBlock-It™ DNA (Rat)	IDBK1004	0.5 mg	Rat repeat sequense DNA, Cytogenetics, Accessories	<input type="checkbox"/>
Mountant and DAPI	CB1172B	50 Test	Antifade & DAPI, Cytogenetics, Accessories	<input type="checkbox"/>
Mouse Metaphase Control Slides (Male Mouse)	AK1001	10 Test	5 Slides/PK, Used for hybridization in FISH, Cytogenetics, Accessories	<input type="checkbox"/>
Sodium Citrate Solution, 0.8%	IS1130	500 ml	Solution 0.8% solution, Cytogenetics, Accessories	<input type="checkbox"/>

Aneuploidy Kits

Mouse IDetect™ Aneuploidy FISH Kit (Chr Y/X/8 cocktail)	IDMFK001	5 Test/Chr	Conc, FISH probe, Mouse Probe, Yred/Xgreen/8yellow(Red+Green)	<input type="checkbox"/>
---	----------	------------	---	--------------------------

Bovine Probe

Bovine IDetect™ Chr X FISH Paint Probe Green	IDBF1061	5 Test	Conc, FISH probe, Cytogenetics, Bovine Probe	<input type="checkbox"/>
Bovine IDetect™ Chr X FISH Probe Red	IDBR1062	5 Test	Conc, FISH probe, Cytogenetics, Bovine Probe	<input type="checkbox"/>
Bovine IDetect™ Chr Y FISH Paint Probe Green	IDBF1058	5 Test	Conc, FISH probe, Cytogenetics, Bovine Probe	<input type="checkbox"/>
Bovine IDetect™ Chr Y FISH Probe Red	IDBR1059	5 Test	Conc, FISH probe, Cytogenetics, Bovine Probe	<input type="checkbox"/>

Product Name**Catalogue****size****Details****NEW**

Bovine IDetect™ PanCentromeric Autosomal FISH Point Probe Biotin	IDBB1071	5 Test	Conc, FISH probe, Cytogenetics, Bovine Probe	<input type="checkbox"/>
Bovine IDetect™ PanCentromeric Autosomal FISH Point Probe Green	IDBF1072	5 Test	Conc, FISH probe, Cytogenetics, Bovine Probe	<input type="checkbox"/>
Bovine IDetect™ PanCentromeric Autosomal FISH Point Probe Red	IDBR1073	5 Test	Conc, FISH probe, Cytogenetics, Bovine Probe	<input type="checkbox"/>

Centromeric Probe

Human IDetect™ Centromeric Chr 1+5+19 Biotin	CB2024	10 T/10 µl	(CONC)	<input type="checkbox"/>
Human IDetect™ Centromeric Chr 1+5+19 Green	CB2025	10 T/10 µl	(CONC)	<input type="checkbox"/>
Human IDetect™ Centromeric Chr 10 biotin	CB1203	10 T/10 µl	(CONC)	<input type="checkbox"/>
Human IDetect™ Centromeric Chr 10 Green	CB1120	10 T/10 µl	(CONC)	<input type="checkbox"/>
Human IDetect™ Centromeric Chr 11 Biotin	CB1122	10 T/10 µl	(CONC)	<input type="checkbox"/>
Human IDetect™ Centromeric Chr 11 Green	CB1121	10 T/10 µl	(CONC)	<input type="checkbox"/>
Human IDetect™ Centromeric Chr 12 Biotin	CB1124	10 T/10 µl	(CONC)	<input type="checkbox"/>
Human IDetect™ Centromeric Chr 12 Green	CB1123	10 T/10 µl	(CONC)	<input type="checkbox"/>
Human IDetect™ Centromeric Chr 13+21 Biotin	CB1204	10 T/10 µl	(CONC)	<input type="checkbox"/>
Human IDetect™ Centromeric Chr 13+21 Green	CB1125	10 T/10 µl	(CONC)	<input type="checkbox"/>
Human IDetect™ Centromeric Chr 14+22 Biotin	CB1127	10 T/10 µl	(CONC)	<input type="checkbox"/>
Human IDetect™ Centromeric Chr 14+22 Green	CB1126	10 T/10 µl	(CONC)	<input type="checkbox"/>
Human IDetect™ Centromeric Chr 15 Biotin	CB1205	10 T/10 µl	(CONC)	<input type="checkbox"/>
Human IDetect™ Centromeric Chr 15 Green	CB1128	10 T/10 µl	(CONC)	<input type="checkbox"/>
Human IDetect™ Centromeric Chr 16 Biotin	CB1207	10 T/10 µl	(CONC)	<input type="checkbox"/>
Human IDetect™ Centromeric Chr 16 Green	CB1206	10 T/10 µl	(CONC)	<input type="checkbox"/>
Human IDetect™ Centromeric Chr 17 Biotin	CB1208	10 T/10 µl	(CONC)	<input type="checkbox"/>
Human IDetect™ Centromeric Chr 17 Green	CB1129	10 T/10 µl	(CONC)	<input type="checkbox"/>
Human IDetect™ Centromeric Chr 18 Biotin	CB1209	10 T/10 µl	(CONC)	<input type="checkbox"/>
Human IDetect™ Centromeric Chr 18 Green	CB1130	10 T/10 µl	(CONC)	<input type="checkbox"/>
Human IDetect™ Centromeric Chr 2 Biotin	CB1109	10 T/10 µl	(CONC)	<input type="checkbox"/>
Human IDetect™ Centromeric Chr 2 Green	CB1108	10 T/10 µl	(CONC)	<input type="checkbox"/>
Human IDetect™ Centromeric Chr 20 Biotin	CB1211	10 T/10 µl	(CONC)	<input type="checkbox"/>
Human IDetect™ Centromeric Chr 20 Green	CB1210	10 T/10 µl	(CONC)	<input type="checkbox"/>

Product Name**Catalogue****size****Details****NEW**

Human IDetect™ Centromeric Chr 3 Biotin	CB1111	10 T/10 µl	(CONC)	<input type="checkbox"/>
Human IDetect™ Centromeric Chr 3 Green	CB1110	10 T/10 µl	(CONC)	<input type="checkbox"/>
Human IDetect™ Centromeric Chr 6 Biotin	CB1202	10 T/10 µl	(CONC)	<input type="checkbox"/>
Human IDetect™ Centromeric Chr 6 Green	CB1113	10 T/10 µl	(CONC)	<input type="checkbox"/>
Human IDetect™ Centromeric Chr 7 Biotin	CB1115	10 T/10 µl	(CONC)	<input type="checkbox"/>
Human IDetect™ Centromeric Chr 7 Green	CB1114	10 T/10 µl	(CONC)	<input type="checkbox"/>
Human IDetect™ Centromeric Chr 8 Biotin	CB1117	10 T/10 µl	(CONC)	<input type="checkbox"/>
Human IDetect™ Centromeric Chr 8 Green	CB1116	10 T/10 µl	(CONC)	<input type="checkbox"/>
Human IDetect™ Centromeric Chr 9 Biotin	CB1119	10 T/10 µl	(CONC)	<input type="checkbox"/>
Human IDetect™ Centromeric Chr 9 Green	CB1118	10 T/10 µl	(CONC)	<input type="checkbox"/>
Human IDetect™ Centromeric Chr X Biotin	CB1134	10 T/10 µl	(CONC)	<input type="checkbox"/>
Human IDetect™ Centromeric Chr X Green	CB1133	10 T/10 µl	(CONC)	<input type="checkbox"/>
Human IDetect™ Centromeric Chr Y Biotin	CB1136	10 T/10 µl	(CONC)	<input type="checkbox"/>
Human IDetect™ Centromeric Chr Y Green	CB1135	10 T/10 µl	(CONC)	<input type="checkbox"/>
Human IDetect™ Pan Centromeric Paint Biotin	CB1576	10 T	(Conc)	<input type="checkbox"/>
Human IDetect™ Pan Centromeric Paint Green	CB1579	10 T	(Conc)	<input type="checkbox"/>

Human Microdeletion Probes

Human IDetect™ DiGeorge N25 FISH Probe	IDR7025	10 Test	N25/RED, Tel Contol/GREEN	<input type="checkbox"/>
Human IDetect™ DiGeorge TUPLE1 FISH Point Probe	IDR7026	10 Test	diGeorge N25/RED, Chr 22 Tel Contol/GREEN	<input type="checkbox"/>

Human WCP

Human IDetect™ Chr 1 Paint Probe Green	IDG7001	10 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 1 Paint Probe Red	IDR7001	10 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 1 Paint Probe Red (Concentrated)	IDR7001-5	5 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 10 Paint Probe Green	IDG7010	10 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 10 Paint Probe Red	IDR7010	10 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 10 Paint Probe Red	IDR7010-5	5 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 11 Paint Probe Green	IDG7011	10 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 11 Paint Probe Red	IDR7011	10 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 11 Paint Probe Red	IDR7011-5	5 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>

Product Name**Catalogue****size****Details****NEW**

Human IDetect™ Chr 12 Paint Probe Red	IDR7012	10 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 12 Paint Probe Red	IDR7012-5	5 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 12 Paint Probe Green	IDG7012	10 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 13 Paint Probe Red	IDR7013	10 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 13 Paint Probe Red	IDR7013-5	5 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 13 Paint Probe Green	IDG7013	10 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 14 Paint Probe Red	IDR7014	10 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 14 Paint Probe Red	IDR7014-5	5 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 14 Paint Probe Green	IDG7014	10 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 15 Paint Probe Red	IDR7015	10 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 15 Paint Probe Red	IDR7015-5	5 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 15 Paint Probe Green	IDG7015	10 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 16 Paint Probe Red	IDR7016	10 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 16 Paint Probe Red	IDR7016-5	5 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 16 Paint Probe Green	IDG7016	10 Test	hWCP, FITC, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 17 Paint Probe Red	IDR7017	10 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 17 Paint Probe Red	IDR7017-5	5 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 17 Paint Probe Green	IDG7017	10 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 18 Paint Probe Red	IDR7018	10 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 18 Paint Probe Red	IDR7018-5	5 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 18 Paint Probe Green	IDG7018	10 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 19 Paint Probe Red	IDR7019	10 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 19 Paint Probe Red	IDR7019-5	5 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 19 Paint Probe Green	IDG7019	10 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 2 Paint Probe Green	IDG7002	10 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 2 Paint Probe Red	IDR7002	10 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 2 Paint Probe Red	IDR7002-5	5 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 20 Paint Probe Red	IDR7020	10 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 20 Paint Probe Red	IDR7020-5	5 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>

Product Name**Catalogue****size****Details****NEW**

Human IDetect™ Chr 20 Paint Probe Green	IDG7020	10 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 21 Paint Probe Red	IDR7021	10 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 21 Paint Probe Red	IDR7021-5	5 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 21 Paint Probe Green	IDG7021	10 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 22 Paint Probe Red	IDR7022	10 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 22 Paint Probe Red	IDR7022-5	5 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 22 Paint Probe Green	IDG7022	10 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 3 Paint Probe Green	IDG7003	10 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 3 Paint Probe Red	IDR7003	10 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 3 Paint Probe Red	IDR7003-5	5 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 4 Paint Probe Green	IDG7004	10 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 4 Paint Probe Red	IDR7004	10 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 4 Paint Probe Red	IDR7004-5	5 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 5 Paint Probe Green	IDG7005	10 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 5 Paint Probe Red	IDR7005	10 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 5 Paint Probe Red	IDR7005-5	5 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 6 Paint Probe Green	IDG7006	10 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 6 Paint Probe Red	IDR7006	10 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 6 Paint Probe Red	IDR7006-5	5 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 7 Paint Probe Green	IDG7007	10 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 7 Paint Probe Red	IDR7007	10 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 7 Paint Probe Red	IDR7007-5	5 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 8 Paint Probe Green	IDG7008	10 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 8 Paint Probe Red	IDR7008	10 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 8 Paint Probe Red	IDR7008-5	5 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 9 Paint Probe Green	IDG7009	10 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 9 Paint Probe Red	IDR7009	10 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr 9 Paint Probe Red	IDR7009-5	5 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr X Paint Probe RED	IDR7023	10 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>

Product Name	Catalogue	size	Details	NEW
Human IDetect™ Chr X Paint Probe RED	IDR7023-5	5 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr X Paint Probe Green	IDG7023	10 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr Y Paint Probe Red	IDR7024	10 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr Y Paint Probe Red	IDR7024-5	5 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Human IDetect™ Chr Y Paint Probe GREEN	IDG7024	10 Test	hWCP, Cytogenetics, Human WCP	<input type="checkbox"/>
Mouse Probe				
Mouse IDetect™ Chr 12 FISH Point Probe (D12Mit197) Biotin	IDMB1052	10 Slide	Conc, FISH probe, Cytogenetics, Mouse Probe	<input type="checkbox"/>
Mouse IDetect™ Chr 19 FISH Point Probe (D19Mit109) Biotin	IDMB1053	10 Slide	Conc, FISH probe, Cytogenetics, Mouse Probe	<input type="checkbox"/>
Mouse IDetect™ Chr 6 FISH Point Probe (D6Mit83) Biotin	IDMB1051	10 Slide	Conc, FISH probe, Cytogenetics, Mouse Probe	<input type="checkbox"/>
Mouse IDetect™ Chr 8 (8.4a) FISH Point Probe Green	IDMF7000	5 Test	Conc, FISH probe for Slide or Suspension, Cytogenetics, Mouse Probe	<input type="checkbox"/>
Mouse IDetect™ Chr 8 (8.5e) FISH Point Probe Green	IDMF7001	5 Test	Conc, FISH probe for Slide or Suspension, Cytogenetics, Mouse Probe	<input type="checkbox"/>
Mouse IDetect™ Chr X FISH Point Probe (DXMit190) Biotin	IDMB1054	10 Slide	Conc, FISH probe, Cytogenetics, Mouse Probe	<input type="checkbox"/>
Mouse IDetect™ Chr X FISH Point Probe (DXMit190) Green	IDMF1074	10 Slide	Conc, FISH probe, Cytogenetics, Mouse Probe	<input type="checkbox"/>
Mouse IDetect™ Chr X FISH Point Probe (DXMit190) Red	IDMR1075	10 Slide	Conc, FISH probe	<input type="checkbox"/>
Mouse IDetect™ Chr X FISH Point Probe Green	IDMF1075	10 Slide	Conc, FISH probe, Cytogenetics, Mouse Probe	<input type="checkbox"/>
Mouse IDetect™ Chr Y FISH Paint Probe Biotin	IDMB1055	10 Slide	Conc, FISH probe, Cytogenetics, Mouse Probe	<input type="checkbox"/>
Mouse IDetect™ Chr Y FISH Paint Probe Green	IDMF1057	10 Slide	Conc, FISH probe, Cytogenetics, Mouse Probe	<input type="checkbox"/>
Mouse IDetect™ Chr Y FISH Paint Probe Red	IDMR1056	10 Slide	Conc, FISH probe, Cytogenetics, Mouse Probe	<input type="checkbox"/>
Mouse IDetect™ Paint Probe Chr 1 Biotin	CB1595M	10 T	(Conc)	<input type="checkbox"/>
Mouse IDetect™ Paint Probe Chr 1 Green	CB1601M	10 T	(Conc)	<input type="checkbox"/>
Mouse IDetect™ Paint Probe Chr 10 Biotin	CB1757M	10 T	(Conc)	<input type="checkbox"/>
Mouse IDetect™ Paint Probe Chr 10 Green	CB1763M	10 T	(Conc)	<input type="checkbox"/>
Mouse IDetect™ Paint Probe Chr 11 Biotin	CB1775M	10 T	(Conc)	<input type="checkbox"/>
Mouse IDetect™ Paint Probe Chr 11 Green	CB1781M	10 T	(Conc)	<input type="checkbox"/>

Product Name**Catalogue****size****Details****NEW**

Mouse IDetect™ Paint Probe Chr 12 Biotin	CB1793M	10 T	(Conc)	<input type="checkbox"/>
Mouse IDetect™ Paint Probe Chr 12 Green	CB1799M	10 T	(Conc)	<input type="checkbox"/>
Mouse IDetect™ Paint Probe Chr 13 Biotin	CB1811M	10 T	(Conc)	<input type="checkbox"/>
Mouse IDetect™ Paint Probe Chr 13 Green	CB1817M	10 T	(Conc)	<input type="checkbox"/>
Mouse IDetect™ Paint Probe Chr 14 Biotin	CB1829M	10 T	(Conc)	<input type="checkbox"/>
Mouse IDetect™ Paint Probe Chr 14 Green	CB1835M	10 T	(Conc)	<input type="checkbox"/>
Mouse IDetect™ Paint Probe Chr 15 Biotin	CB1847M	10 T	(Conc)	<input type="checkbox"/>
Mouse IDetect™ Paint Probe Chr 15 Green	CB1853M	10 T	(Conc)	<input type="checkbox"/>
Mouse IDetect™ Paint Probe Chr 16 Biotin	CB1865M	10 T	(Conc)	<input type="checkbox"/>
Mouse IDetect™ Paint Probe Chr 16 Green	CB1871M	10 T	(Conc)	<input type="checkbox"/>
Mouse IDetect™ Paint Probe Chr 17 Biotin	CB1883M	10 T	(Conc)	<input type="checkbox"/>
Mouse IDetect™ Paint Probe Chr 17 Green	CB1889M	10 T	(Conc)	<input type="checkbox"/>
Mouse IDetect™ Paint Probe Chr 18 Biotin	CB1901M	10 T	(Conc)	<input type="checkbox"/>
Mouse IDetect™ Paint Probe Chr 18 Green	CB1907M	10 T	(Conc)	<input type="checkbox"/>
Mouse IDetect™ Paint Probe Chr 19 Biotin	CB1919M	10 T	(Conc)	<input type="checkbox"/>
Mouse IDetect™ Paint Probe Chr 19 Green	CB1925M	10 T	(Conc)	<input type="checkbox"/>
Mouse IDetect™ Paint Probe Chr 2 Biotin	CB1613M	10 T	(Conc)	<input type="checkbox"/>
Mouse IDetect™ Paint Probe Chr 2 Green	CB1619M	10 T	(Conc)	<input type="checkbox"/>
Mouse IDetect™ Paint Probe Chr 3 Biotin	CB1631M	10 T	(Conc)	<input type="checkbox"/>
Mouse IDetect™ Paint Probe Chr 3 Green	CB1637M	10 T	(Conc)	<input type="checkbox"/>
Mouse IDetect™ Paint Probe Chr 4 Biotin	CB1649M	10 T	(Conc)	<input type="checkbox"/>
Mouse IDetect™ Paint Probe Chr 4 Green	CB1655M	10 T	(Conc)	<input type="checkbox"/>
Mouse IDetect™ Paint Probe Chr 5 Biotin	CB1667M	10 T	(Conc)	<input type="checkbox"/>
Mouse IDetect™ Paint Probe Chr 5 Green	CB1673M	10 T	(Conc)	<input type="checkbox"/>
Mouse IDetect™ Paint Probe Chr 6 Biotin	CB1685M	10 T	(Conc)	<input type="checkbox"/>
Mouse IDetect™ Paint Probe Chr 6 Green	CB1691M	10 T	(Conc)	<input type="checkbox"/>
Mouse IDetect™ Paint Probe Chr 7 Biotin	CB1703M	10 T	(Conc)	<input type="checkbox"/>
Mouse IDetect™ Paint Probe Chr 7 Green	CB1709M	10 T	(Conc)	<input type="checkbox"/>
Mouse IDetect™ Paint Probe Chr 8 Biotin	CB1721M	10 T	(Conc)	<input type="checkbox"/>

Product Name	Catalogue	size	Details	NEW
Mouse IDetect™ Paint Probe Chr 8 Green	CB1727M	10 T	(Conc)	<input type="checkbox"/>
Mouse IDetect™ Paint Probe Chr 9 Biotin	CB1739M	10 T	(Conc)	<input type="checkbox"/>
Mouse IDetect™ Paint Probe Chr 9 Green	CB1745M	10 T	(Conc)	<input type="checkbox"/>
Mouse IDetect™ Paint Probe Chr X Biotin	CB1937M	10 T	(Conc)	<input type="checkbox"/>
Mouse IDetect™ Paint Probe Chr X Green	CB1943M	10 T	(Conc)	<input type="checkbox"/>
Mouse IDetect™ Paint Probe Chr Y Biotin	CB1955M	10 T	(Conc)	<input type="checkbox"/>
Mouse IDetect™ Paint Probe Chr Y Green	CB1961M	10 T	(Conc)	<input type="checkbox"/>
Mouse IDetect™ Pan Centromeric Paint Probe Biotin	CB902M	10 T	(Conc)	<input type="checkbox"/>
Mouse IDetect™ Pan Centromeric Paint Probe Green	CB903M	10 Test	(Conc)	<input type="checkbox"/>
Painting Kit				
Human IDetect™ Chr Paint Kit Biotin	CB1097PD	3 T/Chr	(RTU)	<input type="checkbox"/>
Human IDetect™ Chr Paint Kit Green	CB1098PD	3 T/Chr	(RTU)	<input type="checkbox"/>
Mouse IDetect™ Chr Paint Kit Biotin	CB900M	3 T/Chr	(RTU)	<input type="checkbox"/>
Mouse IDetect™ Chr Paint Kit Green	CB901M	3 T/Chr	(RTU)	<input type="checkbox"/>
Point Probes				
Human IDetect™ Chromosome 13 Point Probe GREEN	IDPPG7013	5 Test	Cytogenetics, Human Point Probes, 13q14.12b-13q14.13b (Chr13 46.1Mb)	<input type="checkbox"/>
Human IDetect™ Chromosome 18 Point Probe AQUA	IDPPA7018	5 Test	Cytogenetics, Human Point Probes, ~18cen (---)	<input type="checkbox"/>
Human IDetect™ Chromosome 21 Point Probe RED	IDPPR7021	5 Test	Cytogenetics, Human Point Probes, 21q22.13b-21q22.2 (Chr21 38.5Mb)	<input type="checkbox"/>
Human IDetect™ Chromosome X Point Probe GREEN	IDPPG7023	5 Test	Cytogenetics, Human Point Probes, Xp11.22-Xp11.2 (ChrX 57.1Mb)	<input type="checkbox"/>
Human IDetect™ Chromosome Y Point Probe RED	IDPPR7024	5 Test	Cytogenetics, Human Point Probes, Yp11.1-Yp11.2 (ChrY 10.7Mb)	<input type="checkbox"/>
Porcine Probe				
Porcine IDetect™ Chr X FISH Paint Probe Green	IDPF1078	5 Test	Conc, FISH probe, Cytogenetics, Porcine Probe	<input type="checkbox"/>
Porcine IDetect™ Chr X FISH Paint Probe Red	IDPR1077	5 Test	Conc, FISH probe, Cytogenetics, Porcine Probe	<input type="checkbox"/>
Porcine IDetect™ Chr Y FISH Paint Probe Green	IDPF1065	5 Test	Conc, FISH probe, Cytogenetics, Porcine Probe	<input type="checkbox"/>
Porcine IDetect™ Chr Y FISH Paint Probe Red	IDPR1065	5 Test	Conc, FISH probe, Cytogenetics, Porcine Probe	<input type="checkbox"/>
Porcine IDetect™ Chr Y FISH Paint Probe Red	IDPR1066	5 Test	Conc, FISH probe, Cytogenetics, Porcine Probe	<input type="checkbox"/>

Product Name**Catalogue****size****Details****NEW****Rat Probe**

Rat IDetect™ Chr X FISH Probe Green	IDRF1067	10 Slide	Conc, FISH probe, Cytogenetics, Rat Probe	<input type="checkbox"/>
Rat IDetect™ Chr X FISH Probe Red	IDRR1068	10 Slide	Conc, FISH probe, Cytogenetics, Rat Probe	<input type="checkbox"/>
Rat IDetect™ Chr Y FISH Paint Probe Green	IDRF1069	10 Slide	Conc, FISH probe, Cytogenetics, Rat Probe	<input type="checkbox"/>
Rat IDetect™ Chr Y FISH Paint Probe Red	IDRR1070	10 Slide	Conc, FISH probe, Cytogenetics, Rat Probe	<input type="checkbox"/>

Telomeric Probe

Human IDetect™ Pan Telomeric Paint Biotin	CB2014	10 T	(Conc)	<input type="checkbox"/>
Human IDetect™ Pan Telomeric Paint Green	CB2015	10 T	(Conc)	<input type="checkbox"/>

Tissue Culture

AmcelGrow® Medium	IS1002	500 ml	To grow amniotic cells (10x500 ml/box), Cytogenetics, Tissue Culture	<input type="checkbox"/>
AmcelGrow® Medium Complete	IS1001	100 ml	To grow amniotic cells, Cytogenetics, Tissue Culture	<input type="checkbox"/>
AmcelGrow® Twin Medium Basal, Less Than 100 bottles	IS1105	90 ml	Must add Supplement cat#IS1106 before use, Cytogenetics, Tissue Culture	<input type="checkbox"/>
AmcelGrow® Twin Medium Basal, Less Than 36 bottles	IS1107	450 ml	Must add Supplement cat#IS1108 before use	<input type="checkbox"/>
AmcelGrow® Twin Medium Supplement, Less Than 100 bottles	IS1106	10 ml	Must use with medium cat#IS1105, Cytogenetics, Tissue Culture	<input type="checkbox"/>
AmcelGrow®Twin Medium Supplement, Less Than 36 bottles	IS1108	50 ml	Must use with medium cat#IS1107	<input checked="" type="checkbox"/>
Amphotericin B Solution, 250 ug / ml	IS1101	100 ml	1-23 Bottles, Similar to Fungizone, Cytogenetics, Tissue Culture	<input checked="" type="checkbox"/>
Amphotericin B Solution, 250 ug / ml, 5mM	IS1066	20 ml	1-23 Bottles, similar to Fungizone (ERSquibb Reg TM), Cytogenetics, Tissue Culture	<input checked="" type="checkbox"/>
Antibiotic / Antimycotic Solution (100X) Liquid	IS1136	100 ml	Contains 10,000 U of penicillin (base), 10,000 µg of streptomycin (base), and 25 µg of amphotericin B/m, in 0.85% saline	<input type="checkbox"/>
Colcemid Solution for Karyotyping	IS1005	10 ml	10µg/ml in PBS for Karyotyping, Cytogenetics, Tissue Culture	<input type="checkbox"/>
Dulbecco's MEM Media 1x Liquid (DMEM)	IS1124-01b	500 ml	1 x concentrate High Glucose Liquid, w L-Glut, 4500 mg /L D-Glucose, 25mM Hepes	<input type="checkbox"/>
Dulbecco's Modified Eagle Medium (DMEM) Powder	IS1018	1 L	Cytogenetics, Tissue Culture	<input type="checkbox"/>

Product Name	Catalogue	size	Details	NEW
Dulbecco's Modified Eagle Medium (DMEM) Powder	IS1019	5 L	(w D-Gluco 4500mg/l, w/o Sod Pyr,w L-Glut, w/o Sod Bicarb), Cytogenetics, Tissue Culture	<input type="checkbox"/>
Dulbecco's Modified Eagle Medium (DMEM) Powder	IS1020	10 L	Cytogenetics, Tissue Culture	<input type="checkbox"/>
Dulbecco's Modified Eagle Medium (DMEM) Powder	IS1021	50 L	Earle's salts base, w L-glut, w/o Sod Bicarb, Cytogenetics, Tissue Culture	<input type="checkbox"/>
Dulbecco's PBS Powder	IS1027	5 L	(w/o Calcium Chl, w/o Magnesium Chl), Cytogenetics, Tissue Culture	<input type="checkbox"/>
Dulbecco's PBS Solution 10 x	IS1124-10	500 ml	Dulbecco's Phosphate Buffered Saline, 10 x concentrate Liquid, w/o Mg, w/o Ca, min 6x500ml	<input type="checkbox"/>
Dulbecco's PBS Solution 1x (D-PBS)	IS1124-01	500 ml	Dulbecco's Phosphate Buffered Saline, 1 x concentrate Liquid, w/o Mg, w/o Ca, min 6x500 ml	<input type="checkbox"/>
Dulbecco's Phosphate Buffered Saline Powder	IS1026	1 L	w/o CaCl ₂ , w/o Mg Cl ₂ , Cytogenetics, Tissue Culture	<input type="checkbox"/>
Dulbecco's Phosphate Buffered Saline Powder	IS1028	10 L	w/o CaCl ₂ , w/o Mg Cl ₂ , Cytogenetics, Tissue Culture	<input type="checkbox"/>
Earle's Balanced Salt Solution (EBSS Complete)	IS1127	500 ml	EBSS complete w Mg, w Ca, min 6x500ml	<input type="checkbox"/>
Gentamycin Sulphate	IS1072	20 ml	1-23 bottles, Tissue Culturte Grade, Cytogenetics, Tissue Culture	<input checked="" type="checkbox"/>
Gentamycin Sulphate	IS1109	100 ml	1-23 Bottles, Liquid, Tissue Culturte Grade, Cytogenetics, Tissue Culture	<input checked="" type="checkbox"/>
Ham's F-10 Nutrient Mixture w L-Glut	IS1114	100 ml	with L-Glutamine, without Sodium Bicarbonate, Cytogenetics, Tissue Culture	<input type="checkbox"/>
Hank's Balanced Salt Solution (HBSS Complete)	IS1126	500 ml	HBSS complete w Mg, w Ca, min 6x500ml	<input type="checkbox"/>
Hanks Balanced Salts (Powder)	IS1118	1 L	Dissolve in 1 Litre of DW, Cytogenetics, Tissue Culture	<input type="checkbox"/>
HEPES Buffer Solution 1M, Lliquid, (1X).	IS1134-1	100 ml		<input type="checkbox"/>
HEPES Buffer Solution 1M, Lliquid, (1X).	IS1134-5	500 ml		<input type="checkbox"/>
KaryoSynKro™ Kit	IS1123	2 x 6 ml (60 T)	4x1.5 ml vials each of 2 solutions, Cytogenetics, Tissue Culture, For use with IS1003/1097/1013	<input type="checkbox"/>
Karyotyping Medium-BM (w/o PHA-M)	IS1012	100 ml	For Bone Marrow Growth, Cytogenetics, Tissue Culture	<input type="checkbox"/>
Karyotyping Medium-BM (w/o PHA-M)	IS1115	500 ml	For Bone Marrow Growth, Cytogenetics, Tissue Culture	<input type="checkbox"/>
Karyotyping Medium-HC (PB&BM) (with conditioned medium)	IS1116	500 ml	For Hematopoietic Cell (PB&BM) Growth, Cytogenetics, Tissue Culture	<input type="checkbox"/>

Product Name	Catalogue	size	Details	NEW
Karyotyping Medium-HC (PB&BM) (with conditioned medium)(100+bottles)	IS1013	100 ml	For Hematopoeitic Cell (PB&BM) Growth, Cytogenetics, Tissue Culture	<input type="checkbox"/>
Karyotyping Medium-PB (w/o PHA-M)	IS1003	100 ml	For Peripheral Blood Cell Growth, Cytogenetics, Tissue Culture	<input type="checkbox"/>
Karyotyping Medium-PB (w/o PHA-M)	IS1011	500 ml	For Peripheral Blood Cell Growth, Cytogenetics, Tissue Culture	<input checked="" type="checkbox"/>
Karyotyping Medium-PB (with PHA-M)	IS1097	100 ml	For Peripheral Blood Cell Growth, Cytogenetics, Tissue Culture	<input type="checkbox"/>
Karyotyping Medium-PB (with PHA-M)	IS1113	5 ml / Tube RTU	RTU TUBE containing 5 ml Peripheral Blood Karyotyping medium with PHA (SAVE cost of using T-25 Flask), Cytogenetics, Tissue Culture	<input type="checkbox"/>
Karyotyping Medium-PB (with PHA-M)	IS1117	500 ml	For Peripheral Blood Cell Growth, Cytogenetics, Tissue Culture	<input type="checkbox"/>
L-Glutamine (200mM) 100x Concentrate in Saline	IS1064	20 ml	1-23 Bottles, Cytogenetics, Tissue Culture	<input checked="" type="checkbox"/>
L-Glutamine (200mM) in Saline Solution	IS1099	100 ml	1-23 Bottles, 100xconc, Cytogenetics, Tissue Culture	<input checked="" type="checkbox"/>
L-Glutamine (200mM) in Saline Solution	IS1100	500 ml	1-23 Bottles, 100xconc, Cytogenetics, Tissue Culture	<input checked="" type="checkbox"/>
MEM non Essential Medium Amino Acids Solution 10mM (100X)	IS1132	100 ml	Liquid, without L-Glutamine. Pack of 100 ml.	<input type="checkbox"/>
MEM Vitamins Solution (100X), Liquid, Pack of 100 ml	IS1133	100 ml		<input type="checkbox"/>
Minimum Essential Medium Eagle (MEM) Powder	IS1014	1 L	Earle's salts base, w L-glut, w/o Sod Bicarb, Cytogenetics, Tissue Culture	<input type="checkbox"/>
Minimum Essential Medium Eagle (MEM) Powder	IS1015	5 L	(Earle's Salts Base, w L-Glut, w/o Sod Bicarb), Cytogenetics, Tissue Culture	<input type="checkbox"/>
Minimum Essential Medium Eagle (MEM) Powder	IS1016	10 L	Earle's salts base, w L-glut, w/o Sod Bicarb, Cytogenetics, Tissue Culture	<input type="checkbox"/>
Minimum Essential Medium Eagle (MEM) Powder	IS1017	50 L	Earle's salts base, w L-glut, w/o Sod Bicarb, Cytogenetics, Tissue Culture	<input type="checkbox"/>
PBS Solution 10x	IS1125-10	500 ml	Phosphate Buffered Saline, 10 x concentrate Liquid, w/o Mg, w/o Ca, min 6x500ml	<input type="checkbox"/>
PBS Solution 1x	IS1125-01	500 ml	Phosphate Buffered Saline, 1 x concentrate Liquid, w/o Mg, w/o Ca, min 6x500ml	<input type="checkbox"/>
Penicillin-Streptomycin Solution 10x	IS1067-10	100 ml	100000U/ml, 100mg/ml in saline, recom diln 1/1000, Cytogenetics, Tissue Culture	<input checked="" type="checkbox"/>
Penicillin-Streptomycin Solution 10x	IS1129-10	20 ml	100000U/ml, 100mg/ml in saline, recom diln 1/1000, Cytogenetics, Tissue Culture	<input checked="" type="checkbox"/>

Product Name	Catalogue	size	Details	NEW
Penicillin-Streptomycin Solution 1x	IS1067-01	100 ml	10000U/ml, 10 mg/ml in saline, recom diln 1/100, Cytogenetics, Tissue Culture	<input checked="" type="checkbox"/>
Penicillin-Streptomycin Solution 1x	IS1128-01	20 ml	10000U/ml, 10 mg/ml in saline, recom diln 1/100, Cytogenetics, Tissue Culture	<input checked="" type="checkbox"/>
Penicillin-Streptomycin-Amphotericin Solution Solution	IS1120	20 ml	1-23 Bottles, 10000U/ml, 10 mg/ml, ??U/ml in saline, recom diln 1/100, Cytogenetics, Tissue Culture	<input checked="" type="checkbox"/>
PHA-M (lyophilized)	IS1007	5 ml	Cytogenetics, Tissue Culture	<input type="checkbox"/>
Potassium Chloride Solution, 0.075M	IS1006	100 ml	(10ug/ml in PBS), ConC: 0.075M, Cytogenetics, Tissue Culture	<input type="checkbox"/>
RPMI 1640-Liquid Media w Hepes	IS1119	100 ml	1-23 Bottles, Contains Hepes 25 mmoles, Cytogenetics, Tissue Culture	<input type="checkbox"/>
RPMI Medium 1640 Liquid	IS1029	500 ml	6+ Bottles, w L-Glut, w Antibiotics, w/o Sod Bicarb, Cytogenetics, Tissue Culture	<input type="checkbox"/>
RPMI Medium 1640 Liquid	IS1098	100 ml	12+ Bottles, w L-Glut, w Antibiotics, w/o Sod Bicarb, Cytogenetics, Tissue Culture	<input type="checkbox"/>
RPMI Medium 1640 Powder	IS1022	1 L	Cytogenetics, Tissue Culture	<input type="checkbox"/>
RPMI Medium 1640 Powder	IS1023	5 L	(w L-Glut, w/o Sod Bicarb), Cytogenetics, Tissue Culture	<input type="checkbox"/>
RPMI Medium 1640 Powder	IS1024	10 L	(w L-Glut, w/o Sod Bicarb), Cytogenetics, Tissue Culture	<input type="checkbox"/>
RPMI Medium 1640 Powder	IS1025	50 L	(w L-Glut, w/o Sod Bicarb), Cytogenetics, Tissue Culture	<input type="checkbox"/>
Sodium Bicarbonate solution 7.5%, Liquid	IS1135	100 ml		<input type="checkbox"/>
Trypsin EDTA 10x (0.5% Trypsin, 0.25% EDTA)	IS1008	20 ml	Cytogenetics, Tissue Culture, dilute1/2 for 0.25% Trypsin , dilute 1/10 for 0.05% Trypsin	<input type="checkbox"/>
Trypsin EDTA 10x in saline (0.5% Trypsin, 0.25% EDTA)	IS1009	100 ml	Cytogenetics, Tissue Culture, dilute1/2 for 0.25% Trypsin , dilute 1/10 for 0.05% Trypsin	<input type="checkbox"/>
Trypsin EDTA Solution A w Phenol Red(0.25% Trypsin,EDTA 0.02%)	IS1137-01-100	100 ml	Cytogenetics, Tissue Culture, Trypsin in Puck's Solution	<input type="checkbox"/>
Trypsin EDTA Solution A w Phenol Red(0.25% Trypsin,EDTA 0.02%)	IS1137-01-500	500 ml	Cytogenetics, Tissue Culture, Trypsin in Puck's Solution	<input type="checkbox"/>
Trypsin EDTA Solution A w Phenol Red(0.25% Trypsin,EDTA 0.05%)	IS1138-01-500	500 ml	Cytogenetics, Tissue Culture, Trypsin in Puck's Solution	<input type="checkbox"/>
Trypsin Solution II (0.25% Trypsin)	IS1008-01	500 ml	Cytogenetics, Tissue Culture, Trypsin	<input type="checkbox"/>
Trypsin Solution II, 10x (2.5% Trypsin)	IS1008-10	100 ml	Cytogenetics, Tissue Culture, dilute 1/10 for 0.25% Trypsin	<input type="checkbox"/>

<p><u>Head Office Address</u> ID Labs Inc. University Research Park 100 Collip Circle, Unit 117 London, Ontario Canada N6G 4X8</p>	<p><u>Mailing Address</u> ID Labs Inc. PO Box 1145 Station CSC London, Ontario Canada N6A 5K2</p>
---	--

Location	Phone	Fax
North America	+1-800-463-4782	+1-800-887-0131
UK / EU	0-800-897-508	0-800-897-509
International	+1-519-434-5057	+1-519-434-2639

www.idlabs.com
[idinfo@idlabs.com](mailto:info@idlabs.com)