


HUMAN IDETECT™ CHROMOSOME FISH POINT PROBES

Chr	Probe	Cytogenetic Location	Label	Cat #	
7	EGFR	Chr 7p12	Aqua	IDPPA7007	NEW!
8	C-MYC	Chr 8q24.21	Green	IDPPG7008	NEW!
9	P16	Chr 9p21	Red	IDPPR7009	NEW!
17	HER-2 / neu	Chr 17q11.2-q12	Orange	IDPPO2017	NEW!
20	ZNF217	Chr 20q13.2	Far Red	IDPPP7020	NEW!

IDetect™ FISH probes in 5 colours from ID Labs™

Colour	Fluorochrome	Ex. (nm)	Em. (nm)
Aqua	IDYE™ 415	418±15	467±10
Green	IDYE™ 495	493±10	521±10
Red	IDYE™ 556	548±15	573±20
Orange	IDYE™ 616	611±10	631±15
Far-red	IDYE™ 647	653±10	672±10



More IDetect™ probes available soon. If you are looking for Custom probes please contact us.

Other Cytogenetics Media and FISH Probes Available

- AmcelGrow® Media for Amniotic Cell Culture
- Karyotyping Media for Hematopoietic, Bone Marrow and Peripheral Blood Cells
- Human IDetect™ Whole Chromosome Paint Probes
- Mouse, Rat, Bovine and Porcine IDetect™ FISH Probes
- Mouse Tri-Color IDetect™ Aneuploidy/Enumeration FISH Probe Kit (Chr 8,X,Y)

Please visit our website at www.idlabs.com.

We can be contacted by email at idinfo@idlabs.com, phone at 519-434-5057 or fax at 519-434-2639.